

Bio

QC Control

P0138

**HBV genotype reference panel
for blood screening**

RUO

REF

P0138



The kit insert contains a detailed protocol and should be read carefully before testing the run control to ensure optimal performance



Table of contents

Intended Use.....	3
Key to Symbols Used.....	3
Summary and Explanation	3
HBV genotypes	3
Principles of the Evaluation Procedure.....	4
Traceability to 1 st WHO standard for HBV-DNA.	4
P0138 HBV 100 copies/ml genotype reference panel for blood screening	5
Storage Instructions	5
Warning and precautions.....	5
Test procedure.....	6
Limitations	6
References	7

www.h-h-c.com

Intended Use

The HBV genotype panel for blood screening tests provides a panel of quantified HBV-DNA preparations covering a large proportion of the currently available genotypes. All members were quantified to obtain dilutions with the same concentration in copies/ml. It can be used to investigate genotype detection efficiency of molecular diagnostic test procedures on Hepatitis B virus DNA in blood samples. This product can be used with amplification methods, including (kinetic) TMA and real-time PCR assays and is useful for development and validation of nucleic acid test systems. It also can be used as a release panel for new batches of HBV-DNA test reagents. This product is not for diagnostic use and for research use only.

Key to Symbols Used



Manufacturer



Lot number



Catalogue number



Store below -30°C



Biological substance
Category B



Research use only



Expiry date



Contents



Caution



Read instructions for use

Summary and Explanation

The HBV genotype panel for blood screening is designed for testing the analytical sensitivity or quantification limits of HBV-DNA tests. The reference panel helps ensure that procedures for HBV-DNA testing are properly validated, and that test results with an unknown group, genotype are consistent across manufacturers, testing laboratories, operators, platforms and assay formats. The bioQControl plasma HBV-DNA virus standards have been diluted in a pool of plasma units that tested negative for viral markers in individual donation NAT and serology testing. All viral standards were quantified by testing in the Siemens Versant bDNA 3.0 assays^{1,2}. The viral concentrations in the plasma pool are ensured by gravimetrically recorded dilutions from calibrated viral stock solutions stored at -70°C.

HBV genotypes

In 1988 Okamoto et al³ divide HBV into 4 genotypes based on a divergence of $\geq 8\%$ in the complete genomic sequence and genotypes A,B,C and D were identified. The relationship between serotypes and genotypes is not clearly known. The same serotype may be classified into different genotypes³. Nordor et al⁴ identified genotypes E and F which differed by more than 4% in the S gene from the other genotypes. Genotype G is reported in 2000 from samples of French and American patients⁵ but its geographic origin is still unknown⁶. The precore and core regions of genotype G are aberrant with a 36-nucleotide insertion within the core gene making it the longest of the HBV genotypes⁷. Genotype I

described in Vietnam⁸ may not meet the criteria for a novel genotype since the diversity in its complete genome sequence is only 7% from that of its closest neighbour, genotype C⁹. Genotype J is a novel variant described in a Japanese patient. It is thought to be phylogenetically positioned between human and primate HBV variants being close to strains which had been previously found in orang-utans and gibbons¹⁰. Subgenotypes are also described if there is a divergence of > 4% (but less than 7.5%) of the nucleotide sequence in the complete genomic sequence. Divergence of < 4% between subgenotypes are referred to as "clades".

Principles of the Evaluation Procedure.

The HBV genotype panel for blood screening members have been carefully formulated to mimic human plasma specimens containing 100 copies/ml (19 IU/ml) HBV-DNA. The HBV genotype panel for blood screening is suitable for evaluate the ability of the assay specific primers and probes to recognise all HBV genotypes. The composition of the panel covers world-wide most spread HBV-variants. As HBV is continuously evolving we recognise not all variants are included. Laboratories should find equal, positive responses for the different samples.

Traceability to 1st WHO standard for HBV-DNA.

A standard dilution of HBV-DNA genotype A preparation was included in the WHO collaborative study¹¹ to establish the 1st and 2nd WHO standard for HBV-DNA. Later we calibrated the genotype A standard against the 1st WHO standard by testing standard dilution series in bDNA 3.0. It was found one bDNA copy is equal to 5.33 IU of the first WHO HBV-DNA standard. The other genotypes were calibrated on the genotype A preparation (sample 01). All panel members contain 19 IU/ml.

P0138 HBV 100 copies/ml genotype reference panel for blood screening

Each panel member is quantified at 100 copies/ml[#] and filled off with 4.3 ml.

Table 1 composition of the panel

Member	HBV genotype	HBV serotype	Country of origin
1	A	adw2	Netherlands
2	B	ayw1	Indonesia
3	C	adr	USA
4	D	ayw2	USA
5	E	ayw3	USA
6	F	adw4	USA
7	G	adw2	USA
8	A2	adw2	Germany
9	D	ayw2/3	Germany
10	A1	adw2	South Africa
11	A1	adw2	Brasilia
12	A2	adw2	Germany
13	B1	adw2	Japan
14	B2	adw2	Japan
15	B4	ayw1	Vietnam
16	C2	adr	Japan
17	C2	adr	Japan
18	C2	adr	Russia
19	D1	ayw2	Germany
20	D3	ayw2	South Africa
21	D1	ayw3	Iran
22	E	ayw4	West Africa
23	F3	adw4	Brasilia
24	G	adw2	Germany
25	negative		

The HBV-DNA standards has been diluted in a pool of plasma units that tested individually negative for HBsAg, anti-HBc, anti-HBs, anti-HCV, anti-HIV 1 and 2, HBV-DNA, HCV-RNA and HIV-1 RNA.

Storage Instructions

It is recommended to store the panel at -30°C or lower to ensure highest quality. Discard any unused material after the first use. Any panel members that appear cloudy or contain precipitates after thawing should be discarded.

Warning and precautions

The P0138 HBV 100 copies/ml genotype reference panel members contain infectious *HBV* virions and are bio-hazardous. Observe the universal precautions for prevention of transmission of infectious agents when handling these materials^{12,13,14}. Although the normal human plasma used in the production of this panel was negative for infectious disease markers the reference panel members should be handled as if capable of transmitting (unknown) infectious agents.

- Do not pipette by mouth.
- Use personal protective equipment, including lab coats, gloves and safety glasses.
- Do not eat, drink or smoke in areas where P0138 HBV 100 copies/ml genotype reference panel and specimens are handled.
- Disinfect liquids, materials or spills with a 0.5% sodium hypochlorite solution or equivalent.
- Dispose of all materials and liquids used in the procedure as if they contained pathogenic agents.

Test procedure

- Thaw the panel members quickly in a water bath at 37°C to avoid formation of cryo-precipitates.
- Mix gently during thawing until ice clot has disappeared.
- Immediately after thawing, vortex briefly, and give a short spin before releasing screw cap from the vials.
- The panel members should be handled and tested in a manner identical to that required for clinical specimens run in the test procedure being evaluated.
- Follow the manufacturers or testing laboratory instructions and recommendations for the handling and testing of clinical specimens.

Limitations

The P0138 HBV 100 copies/ml genotype reference panel is not intended to replace the internal calibrators integral to in vitro diagnostic (IVD) test kits, but may be used as external, independent standards for the assessment of the performance of qualitative or quantitative NAT assays. The panel is not an in vitro diagnostic and for research use only.

References

1. Package insert Siemens Versant HBV-DNA bDNA 3.0 assay
2. Collins ML, Zayati C, Detmer JJ, Daly B, Kolberg JA, Cha TA, Irvine BD, Tucker J, Urdea MS Preparation and characterization of DNA standards for use in quantitative branched DNA hybridization assays. *Anal Biochem.* 1995 Mar 20;226(1):120-9
3. Okamoto H, Tsuda F, Sakugawa H, Sastrosoewignjo RI, Imai M, Miyakawa Y, Mayumi M. Typing hepatitis B virus by homology in nucleotide sequence: comparison of surface antigen subtypes. *J Gen Virol.* 1988;69:(Pt 10):2575-2583.
4. Norder H, Hammas B, Löfdahl S, Couroucé AM, Magnius LO. Comparison of the amino acid sequences of nine different serotypes of hepatitis B surface antigen and genomic classification of the corresponding hepatitis B virus strains. *J Gen Virol.* 1992;73:(Pt 5):1201-1208
5. Stuyver L, De Gendt S, Van Geyt C, Zoulim F, Fried M, Schinazi RF, Rossau R. A new genotype of hepatitis B virus: complete genome and phylogenetic relatedness. *J Gen Virol.* 2000;81:67-74.
6. Lindh M. HBV genotype G-an odd genotype of unknown origin. *J Clin Virol.* 2005;34:315-316.
7. Li K, Zoulim F, Pichoud C, Kwei K, Villet S, Wands J, Li J, Tong S. Critical role of the 36-nucleotide insertion in hepatitis B virus genotype G in core protein expression, genome replication, and virion secretion. *J Virol.* 2007;81:9202-9215.
8. Tran TT, Trinh TN, Abe K. New complex recombinant genotype of hepatitis B virus identified in Vietnam. *J Virol.* 2008;82:5657-5663.
9. Kurbanov F, Tanaka Y, Kramvis A, Simmonds P, Mizokami M. When should "I" consider a new hepatitis B virus genotype?. *J Virol.* 2008;82:8241-8242.
10. Tatematsu K, Tanaka Y, Kurbanov F, Sugauchi F, Mano S, Maeshiro T, Nakayoshi T, Wakuta M, Miyakawa Y, Mizokami M. A genetic variant of hepatitis B virus divergent from known human and ape genotypes isolated from a Japanese patient and provisionally assigned to new genotype J. *J Virol.* 2009;83:10538-10547.
11. An international Collaborative study to establish a WHO International Standard for Hepatitis B virus DNA Nucleic Amplification Techniques. J. Saldanha, W.H. Gerhlich, P.N. Lelie, K.H. Heerman and A. Heath. *Vox Sang.* 2001, 80: 63-7
12. Centers for Disease Control (CDC). Recommendations for prevention of HBV transmission in health care settings. *MMWR* 1987; 36 (supplement no. 2S).
13. Centers for Disease Control (CDC). Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in health-care settings. *MMWR* 1988; 37:377-388.
14. Centers for Disease Control (CDC). Guidelines for prevention of transmission of human immunodeficiency virus and hepatitis B virus to health-care and public-safety workers. *MMWR* 1989; 38(S-6): 1-36.

www.h-h-c.com



BioQControl B.V.
Visseringlaan25
2288 ER Rijswijk
The Netherlands

Tel: +31 (0)88 235 33 33
Fax: +31 (0)88 235 33 00
Internet: www.BioQControl.com

KI4138
V1.0 November 2015