

PRODUCT DESCRIPTION:

Severe acute respiratory syndrome-related coronavirus 2 (SARS-CoV-2), isolate hCoV-19/USA/COR-22-063113/2022 was isolated from a human in May 2022, in Memphis, Tennessee, USA. Under the nomenclature system introduced by GISAID (Global Initiative on Sharing All Influenza Data), SARS-CoV-2, isolate hCoV-19/USA/COR-22-063113/2022 is assigned lineage BA.5 using the Phylogenetic Assignment of Named Global Outbreak lineages (PANGO) tool. It was labeled as an Omicron variant by the World Health Organization (WHO).

SARS-CoV-2 was propagated in VeroE6 cells. The viral culture fluid was clarified via centrifugation to remove cellular debris prior to inactivation.

PROtrol™ (Protein control) is formulated with viral particles that have been modified to render them non-infectious. Each vial contains 1.0 mL of SARS-CoV-2 PROtrol™ frozen in 2% MEM media.

The pre-inactivation titer was determined from an infectious aliquot. Each lot of PROtrol™ is tested using an N-protein ELISA developed internally at ZeptoMetrix.

Viral inactivation is verified by the absence of viral growth in tissue culture-based infectivity assay.

INTENDED USE:

PROtrol™ is intended for use in analytical and quality control testing of antigen-based assays. The suitability and performance characteristics should be determined by your laboratory for each intended usage.

These products are NOT intended for use in the manufacture or processing of injectable products subject to licensure under section 351 of the Public Health Service Act or for any other product intended for administration to humans.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

BIOSAFETY:

Please consult your institution's regulations regarding the use of this product. For a detailed discussion on biological safety see the current edition of Biosafety in Microbiological and Biomedical Laboratories (BMBL), published by the CDC.

PRECAUTIONS:








- Use Universal Precautions, this product is potentially biohazardous.
- Repetitive freezing and thawing is not recommended.
- To avoid cross-contamination, use separate pipette tips for all reagents.

RECOMMENDED STORAGE:

PROtrol™ should be stored at -65°C or below.

¹The following reagent was obtained through BEI Resources, NIAID, NIH: SARS-Related Coronavirus 2, Isolate hCoV-19/USA/COR-22-063113/2022 (Lineage BA.5; Omicron Variant), NR-58616, contributed by Dr. Richard J. Webby.

PIPROSARS(COV2)-658
Revision: 00
Effective Date: 03/20/2023

 REF	Catalog Number		Temperature Limitation
 LOT	Batch Code		Expiration Date
 RUO	For Research Use Only		Biological Risk
	Manufacturer		

www.ZeptoMetrix.com

ZeptoMetrix LLC • 878 Main Street, Buffalo, NY 14202 USA • Tel (800) 274-5487
This product was manufactured in a facility which has a Quality Management System that is ISO 13485 certified.

PCA# N/A
Page 1 of 1